

ALUMINIUM BOX - 17.5X7.6X2 cm



Product Code	26662
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2740347
CND	V0499
EAN/UPC	8023279266627
GTIN	08023279266627
GMDN	13730

Description

ALUMINIUM BOX 17.5x7.6xh 2 cm

Useful containers for the sterilization of surgical instruments and to maintain them sterile.
Silicone tray not included.
Nice aluminium box with blue base and silver cover.