



Product Code	36666
Unit of sale	box of 600 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	802327936662
PARAF	977815897

Description

BIBS WITH POCKET

Disposable bibs made of 1 ply of highly absorbent wadding of cellulose and 1 ply of waterproof blue polyethylene.

Box of 600 (6 bags of 100).

Size: 37 x 67 cm + 10 cm (pocket)
 Easy fastening on patient by lacets
 Extremely soft and comfortable
 Large pocket