

MESOTHERAPY LUER NEEDLES 31G 0,26x4 mm - light blue



Product Code	23674
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2648224
CND	A019006
EAN/UPC	8023279236743
PARAF	951615436

Description

MESOTHERAPY LUER NEEDLES - STERILE

- Gauge 31G (Ø 0.26 mm)
- Length: 4 mm
- Colour: light blue
- Box of 100 pcs.
- Sterile

Meso-relle disposable hypodermic needles are multiple sharpened to facilitate penetration. Stainless steel AISI 304 cannula and polypropilene connection for syringe.

Sterilized with ethylene oxide.
Sold in individual blister made of medical paper.

Made in Japan.