

HEYMANN SCISSORS - 18 cm



| | |
|---------------|----------------|
| Product Code | 26781 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | No |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 142779 |
| CND | L01040602 |
| EAN/UPC | 8023279267815 |
| PARAF | 903600373 |
| GTIN | 08023279267815 |
| GMDN | 38727 |

Description

Stainless Steel Heymann scissors - 18 cm