

ROSEN EAR SUCTION CANNULA diameter 1.5 mm



Product Code	26797
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2847276
CND	L030101
EAN/UPC	8023279267976
PARAF	903149639
GTIN	08023279267976
GMDN	36269

Description

Rosen ear suction cannula Ø 1.5 mm (15.5 cm) with adaptor with LUER tip and thumb valve