

MARTIN HYSTEROMETER - malleable - 30 cm



| | |
|---------------|----------------|
| Product Code | 26816 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | No |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2829392 |
| CND | L99 |
| EAN/UPC | 8023279268164 |
| GTIN | 08023279268164 |
| GMDN | 13644 |

Description

Martin hysteroscope - 30 cm - malleable

Brass handle, copper wire.