

# PARAGON DISPOSABLE SCALPEL BLADES N.20 - sterile



Product Code	26905
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2272184
CND	V010302
EAN/UPC	5033955103064
PARAF	900115559

## Description

"PARAGON" DISPOSABLE S/S BLADES - sterile

- Blades 20

High quality blades manufactured with best British stainless steel guarantee excellent sharpening and perfect cut control.  
Individually packed in metal foil.  
Gamma rays sterilized.

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE

Made in the UK.