


DIAMANTINE STAINLESS STEEL BLADES N. 10 - sterile



Product Code	27010
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2391015
CND	V010101
EAN/UPC	8023279270105

Description

"GIMA" DISPOSABLE BLADES - sterile - dispenser box of 100

- Diamantine Stainless steel

DIAMANTINE™ EDGE HIGH PERFORMANCE BLADES

High quality, super sharp blades with zirconium nitride coating fused on the cutting edge to satisfy the most demanding surgeons requirements.

Durable cutting edge, that stays sharp incision after incision.

The diamond like material creates a very hard edge that will not dull even under the most demanding cutting needs.

Bio-compatible and highly corrosion resistant.

Compatible with autoclavable scalpel handles N. 3 and N. 4.

Stainless steel, sterile.

Manufactured in accordance with requirements of British Standards BS 2982:1992, ISO 7740:1985, ISO 7153-1: 1991, EN 10088-3: 2005.

Each blade is hallmarked and individually packed in aluminium foil.

Blade size is clearly identified on the reverse side of the foil peel-pack and on product unit boxes.

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, HU, SK, SE, FI, DK, RO, BG, GR, HR, SI, Arabic.

MDR Certified.