

CONTAINER WITH FILTER small h 200 mm - blue - perforated



Product Code	37033
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 8023279370331

Description

Standard container with two filters and perforated bottom.

Perforated bottom allows much quicker sterilization process (half the time compared to not perforated). Lock nail on bottom side made of chrome.

Standard containers (supplied with 1 paper filter and filter holder)

- Small model: 285 x 280 mm
- Height: 200 mm
- Perforated bottom
- 2 filters
- Colour: blue

All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.