

CONTAINER WITH FILTER medium h 100 mm - grey



|               |                |
|---------------|----------------|
| Product Code  | 37050          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |
|               |                |
| EAN/UPC       | 8023279370508  |

Description

Standard container with one filter and non perforated bottom.

- Medium model: 465 x 280 mm
- Height: 100 mm
- Non perforated bottom
- 1 filter
- Colour: grey

Standard containers (supplied with 1 paper filter and filter holder).  
All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.