



Product Code	27055
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2247663
CND	V010302
EAN/UPC	0707387789008
PARAF	951090036

Description

3M™ PIVOTING BLADES - 9660 - for Surgical clipper (code 27054)

Pivoting head of 3M™ surgical clipper can works with this kind of blades.

Single use clipper blade individually sealed.

Technical Specifications

- Features of the clipper (code 27054) with PIVOTING blades 9660
- Pivoting head
 - Ergonomically designed handle that's easy to hold
 - Easy-to-attach blade
 - Cuts through all types of head and body hair, wet or dry, in a single pass
 - Away-from-the-skin cutting blade minimizes potential for nicks or cuts, which helps reduce the risk of SSIs (surgical site infections)
 - Li-ion battery designed to hold a charger longer and battery indicator light
 - Approximately 160 minutes of run time on a four-hour recharge
 - Submersible for easy cleaning (meets IP X7 rating of submersible in water no deeper than 1 m for up to 30 min)