



Product Code	27087
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	4015672104945
PARAF	987476518

Description

OMRON HN286 DIGITAL SCALE

Electronic bathroom scale with four integrated sensors for fast and reliable reading, up to 180 kg.
Automatic shutdown after 6 seconds.
Switch on by contact on the tray.

Multilingual manual in all European languages and NO, TR, RU, Arabic.

Technical Specifications

Capacity: 50 to 180 kg
Sensitivity: 100 g
Power: CR2032 lithium battery - included
Size: 300x300xh 18 mm
Weight: 1.9 kg