

# CONTAINER WITH FILTER large h 150 mm - red



Product Code	37112
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2811854
CND	S9099
EAN/UPC	8023279371123

## Description

Standard container with two filters and non perforated bottom.

- Large model: 580 x 280 mm
- Height: 150 mm
- Non perforated bottom
- 2 filters
- Colour: red

Standard containers (supplied with 1 paper filter and filter holder).  
All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.