

"SUPERTEXAS LEATHER" MEDICAL BAG - chestnut



| | |
|---------------|-------------------|
| Product Code | 27116 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | Not required |
| CH-REP | Not required |
| EAN/UPC | 8023279271164 |

Description

With outer pocket and inside zip fastener pocket, 2 burnished locks, inside transversal partition 8 cm from side, large central area.

External dimensions: 41x20xh 22 cm

Made in Italy.