

# "UTILITY" DOCTOR BAG - nylon blue



|               |                   |
|---------------|-------------------|
| Product Code  | 27125             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 8023279271256 |

## Description

### UTILITY DOCTOR BAG - blue - empty

Practical doctor bag for easy transport of instruments and medications, provided with 8 internal movable divisions, 2 removable external pockets and 1 removable mini bag (17x3xh 24 cm), 1 fixed external pocket with zip including 4 elastics retaining and other small compartments.

Size: 36x22xh 25 cm.

Made of water resistant 600D polyester  
Removable minibag