

"UTILITY" DOCTOR BAG - nylon beige

*accessori non inclusi
accessories not included
accessoires non inclus*



Product Code	27127
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279271270
---------	---------------

PARAF	903604116
-------	-----------

Description

UTILITY DOCTOR BAG - beige - empty

Practical doctor bag for easy transport of instruments and medications, provided with 8 internal movable divisions, 2 removable external pockets and 1 removable mini bag (17x3xh 24 cm), 1 fixed external pocket with zip including 4 elastics retaining and other small compartments.

Size: 36x22xh 25 cm.

Made of water resistant 600D polyester
Removable minibag