

PROFESSIONAL BAG

accessori non inclusi
accessories not included
accessoires non inclus



Product Code 27128

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279271287

Description

PROFESSIONAL BAG - empty

Versatile bag, designed for general practitioners and any professional. It includes: 3 detachable transparent pouches with velcro, 1 ampoule holder with 23 elastic bands with velcro, 1 large pouch with 4 dividers with velcro, clip paper board, document holder, various pockets. Shoulder strap, handle, hidden straps to carry on back.

Made in 70D nylon (external) 600D polyester (internal).

Weight: 3.3 kg.
Maximum load: 5 kg
External size: 42x35x15 cm.