

# CONTAINER WITH FILTER large h 135 mm - blue - perforated



Product Code	37131
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2811854
CND	S9099
EAN/UPC	8023279371314

## Description

Standard container with 4 filters and perforated bottom.

Perforated bottom allows much quicker sterilization process (half the time compared to not perforated). Lock nail on bottom side made of chrome.

Standard containers (supplied with 1 paper filter and filter holder)

- Large model: 580 x 280 mm
- Height: 135 mm
- Perforated bottom
- 4 filters
- Colour: blue

All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.