

MULTIUSE SOFT CASE - black



Product Code	27137
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279271379
PARAF	987469196

Description

MULTIUSE SOFT CASE - black

Inner insert divides this bag into two compartments, one for stethoscope, one for nurse accessories. Bottom compartment with one elastic pocket and one transparent zipper pocket. Upper compartment with 6 elastic loops, one transparent zipper pocket and one mesh pocket. Double zip for easy access.
Size: 28x12.5x6.5 cm.
Supplied with carrying handle.