

TRANSPORT BAG - nylon black for ECG, Diatermo, ...



Product Code	27139
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279271393
PARAF	903141137

Description

DIATERMO BAG - black - 600D polyester

Designed to carry small electrosurgical units, audiometers, spirometers, ECG with their accessories.
One removable partition to separate device from accessories and one zip pocket.

Size: 30 x 40 x h14 cm