

# SQUEEZE EGG - X-soft - yellow



|               |                   |
|---------------|-------------------|
| Product Code  | 47140             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 8023279471403 |
| PARAF   | 973923307     |

## Description

### SQUEEZE EGG - X-soft

These TPE (Thermoplastic elastomers) eggs offer variable resistance training for hands and fingers. Help grip strengthening and reduce stress.