

POCKET ORGANIZER - lilac



Product Code	27143
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279271430
---------	---------------

Description

Pocket organizer - lilac

Pocket size organizers help professionals to have easy and fast access to needed tools. Multiple pockets can hold pens, thermometer, penlights.

Size: 15x12xh 2 cm