

SMART BAG - medium - blue



Product Code	27152
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279271522
PARAF	977811797

Description

Smart bag in polyester PVC coated - empty

- Material: polyester PVC coated
- Model: medium
- Colour: blue
- Size: 55 x 35 x h 32 cm
- Internal modular spaces: 6
- External pockets: 4

Professional emergency bags extremely capacious with internal adjustable compartments which permit personalized organization of the space.
Provided with external pockets, internal divisions, and one removable pouch with transparent window.
Bags are externally made of strong, water resistant 600D polyester PVC coated.
All bags are internally made of 600D polyester PVC coated.
Provided with double yellow reflective strips.
Internal lining is removable for easy cleaning.
Waterproof rubber bottom to prevent water infiltration.