

SMART BAG PVC coated - medium - black



Product Code	27158
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279271584
PARAF	977811835

Description

Smart bag in PVC - empty

- Polyester PVC tarpaulin
- Model: medium
- Colour: black
- Size: 55 x 35 x h 38 cm
- Internal modular spaces: 6
- External pockets: 5

Professional emergency bags extremely capacious with internal adjustable compartments which permit personalized organization of the space.
Provided with external pockets, internal divisions, and one removable pouches with transparent window.
Bags are externally made of cold resistant PVC tarpaulin.
All bags are internally made of 600D polyester PVC coated.
Provided with double yellow reflective strips.
Internal lining is removable for easy cleaning.
Waterproof rubber bottom to prevent water infiltration.