

XION FLEXIBLE NOSEPHARINGOSCOPE 18.000 pixel - without cable



Product Code	71606
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2548757
CND	Z12021006
EAN/UPC	8023279716061

Description

XION FLEXIBLE NOSEPHARYNGOSCOPE WITHOUT CABLE

Nosepharyngoscope ideal for endoscopic examination of the nasopharynx, nasal and oral cavities, as well as the trachea.

The high precision multi-lens optics creates an optimal image.

The thin insertion tubes are exceptionally maneuverable and have a special wide angle lens whereby a precise tip control and a remarkable depth of focus are reached.

High quality materials guarantee excellent and homogeneous illumination.

The handpiece is ergonomically formed and the insertion tube is especially comfortable for the patient.

Wide angle, multi-lens objectives for excellent image quality.

Excellent ergonomics allows one hand operation.

Suitable for mechanical cleaning and disinfection.

For illumination, it needs light source (30803), adapter (30830), F.O. cable (see page 239), adapter (30840).

Supplied in lightweight rigid case with leakage tester and user manual (GB, IT). FR, DE manual available on www.gimaitaly.com

Made in Germany.

Technical Specifications

RESOLUTION: 18,000 pixels
WORKING LENGTH: 320 mm
WORKING Ø: 3.4 mm
TIP DEFLECTION: up 130°, down 130°
DIRECTION OF VIEW: 0°
DEPTH OF FOCUS: 1-50 mm
WEIGHT: 230 g
FIELD OF VIEW: 80°