

LIFE-2 BAG



|               |                   |
|---------------|-------------------|
| Product Code  | 27161             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|        |              |
|--------|--------------|
| UK-REP | Not required |
| CH-REP | Not required |

|         |               |
|---------|---------------|
| EAN/UPC | 8023279271614 |
| PARAF   | 977811847     |

Description

LIFE-2 EMERGENCY BAG with 5 coloured pouches

Two section first-aid bag with double opening zip: bottom side with shockproof inner basin, upper side with 5 removable coloured pouches (model E1) with transparent window.

Made of strong, water resistant 600D polyester material and provided with reflective yellow strips.  
Shoulder or hand portable.

Size: 49 x 28 x h 33 cm.