

CUBO BAG PVC coated - red



Product Code	27163
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279271638
PARAF	977811850

Description

CUBO EMERGENCY BAG - PVC tarpaulin - red

Portable emergency bag with internal and external zip pocket and 2 removable pouches with transparent window.

Interior made of strong, water resistant 600D polyester with outer PVC.
Provided with double reflective yellow strip.

Size: 28 x 34 x 13 cm.