

EMERGENCY BAG PVC coated - blue



Product Code	27167
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279271676
PARAF	977811886

Description

EMERGENCY BAG - PVC - blue - 35x45xh 21 cm
Capacious bag, despite the compact size. Large central compartment, 1 front and 3 lateral pockets.
May be carried on shoulder (braces) or hand (2 handles).
Made of strong, water resistant material 600D polyester PVC tarpaulin.
Provided with yellow reflective strips and rubber bottom.

2 lateral pockets 22x18 cm and 11x18 cm
Lateral pocket 32x18 cm
Front pocket 35x20 cm