

LOGIC 2 RUCKSACK



Product Code	27171
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279271713
PARAF	977811924

Description

LOGIC-2 RUCKSACK - polyester - large - 40x25xh 47 cm

Rucksack externally made of strong water resistant 600D polyester PVC coated.
Supplied with 6 coloured pouches with transparent window and velcro strips, 2 frontal pockets and 2 lateral pockets.
Internal velcro strips and elastics, to retain instruments.
Central division with zip and elastics.
Provided with yellow reflective strips and shoulder straps.