

VINYL GLOVES - POWDER FREE ELASTICISED - large



Product Code	25718
Unit of sale	box of 100 pcs
Minimum order	1
Type	PPE
Class	III

UK-REP	Not required
EAN/UPC	8016629000798
PARAF	933882843

Description

VINYL GLOVES - POWDER FREE

- Elasticized for extra sensitivity and optimal grip
- Large size
- NON sterile

Gloves made of polyvinyl chloride, latex-free.
Powder free, ideal for sensitive skins.
For everyday use, they ensure protection and sensitivity.
Suitable for dressing in hospitals and laboratory use.

Packed in dispenser box of 100 units.
MD class I.
PPE cat II.
EN455-EN374-AQL<1.5