

# E6 POUCH with window for Smart Pediatric Bag - grey



|               |                   |
|---------------|-------------------|
| Product Code  | 27200             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 8023279272000 |

## Description

Pouche provided with a large strip of female velcro on back.  
Made in 600D polyester.

Model: E6  
Size: 38 x 10 x 5 cm  
Colour: gray