

THERMAL BAG - yellow nylon



Product Code	27209
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279272093
PARAF	982516534

Description

THERMAL BAG - yellow - nylon

Crafted, durable and insulated thermal bag with detachable shoulder strap.
Supplied with pockets and 3 temperature packs, which maintain cold and hot temperature for long period.

External size: 34x25x24 cm - Internal size: 30x23x21 cm

Made in Italy.