

B.S. CONTAINER WITH VALVE medium h 100 mm - red



Product Code	37210
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC	8023279372106
---------	---------------

Description

Barrier System container with valve and non perforated bottom.
This system is preferred when is requested continuous use of containers, it avoids filter replacement. Valve opens during sterilization to allow exit of steam.
Lock nail on bottom side made of chrome.

- Medium model: 465 x 280 mm
- Height: 100 mm
- Non perforated bottom
- Barrier System with valve
- Colour: red

All containers are autoclavable at maximum 134°C and 2.2 bar.

NOT FOR SALE IN GERMANY.