



Product Code	27233
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	
UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279272338
PARAF	950834224

Description

SPORTPIX SMART WATCH

Smart watches designed to monitor different health and sport parameters.

Sportpix is more suitable for sport activities.

11 PARAMETERS: BODY TEMPERATURE, HEART RATE, ECG, BPM, SpO2, RESPIRATION RATE, SLEEP, SPORTS, CALORIES, STEPS, WEATHER.

In addition to the above parameters, both watches include stopwatch, music player and photos.
Wireless connection to SIMPLO App, compatible with Android 8.1 and above, iOS 13 and above.

Supplied with power cable and manual in GB, IT. FR, ES, DE available on www.gimaitaly.com.

Technical Specifications

Accelerometer: 3D sensor
Display: 1.3" OLED
Data sync: 4.0 BT
Temperature range: -10° ~ 50°
Battery: Lithium polymer
Battery capacity: 240 MAH
Battery life stand-by: 15 days
normal use: 5 days
frequent use: 3 days
Charging time: 2-3 hours
Weight: 38 g
Water resistance: IP67