

GIMA BIG DIAL ANALOGUE SCALE



Product Code	27239
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279272390
PARAF	962663643

Description

GIMA BIG DIAL SCALE

Mechanical scale with classic design: a large and clear dial, a robust housing made of steel.

Supplied in multilingual gift box.

Made in Hungary.

Technical Specifications

Capacity: 160 kg
Sensitivity: 1,000 g
Minimum measuring weight: 5 kg
Dial Ø: 180 mm
Weight: 3 kg
Dimension: 430x380x70 mm