

SECA 750 SCALE



Product Code	27241
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No
RDM (NSIS)	1403012
CND	Y0399
EAN/UPC	4012030012683
PARAF	987475858

Description

Mechanical flat scale with sturdy, powder coated steel casing.
Pretty and stable with its white powder coated lacquer.
A hard wearing synthetic white leather covers weighing platform.

Made in Europe.
User manual in: GB, FR, IT, ES, DE, PT, PL.

Technical Specifications

Capacity: 150 kg - 330 lbs
Sensitivity: 1 kg
Size: 317x117x477 mm
Weight: 3.5 kg