

SOEHNLE 7808 DIGITAL MULTIFUNCTION SCALE



Product Code	27243
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	4006484774198

Description

MULTIFUNCTIONAL SCALE - SOEHNLE 7808
Manual in: DE, GB, FR, NL, ES, IT, GR, FI, NO, SE, DK, CZ, RU.

Multifunction scale prepared for easy (and later) mounting of support handrails, adiposity handrails, bracket for indicator, drive-on ramps.
Movable by means of 2 rollers and handle.

Large LCD display, digit height 20 mm.

Made in Germany.

Technical Specifications

Capacity: 300 kg
Sensitivity: 100 g
Dimension: 840x880x55 mm
Weight: 16 kg
Power supply: rechargeable battery with charger (on request power supply 230V - 50-60 Hz)