

GLASS DIGITAL SCALE - black



Product Code	27248
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279272482
PARAF	976326987

Description

DIGITAL GLASS SCALE - black

Electronic bathroom scale with 6 mm thick tempered glass platform.
Ultra slim design: only 19 mm height.
4 high precision G sensors.
Automatic on and off function.
Lb/kg/st switchable by user.

LCD display 72x35 mm

Multilingual colour box and manual: GB, FR, IT, ES, DE, PL, PT, HU, RO, FI, SE, HR, GR, Arabic

Technical Specifications

Capacity: max 150 kg (330 lb) / min. 2 kg (4.4 lb)
Sensitivity: 100 g (0.2 lb)
Dimension: 300x300x19 mm
Power: 1xCR2032 lithium battery - not included
Warranty: 2 years