



Product Code	27249
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279272499
PARAF	976326975

## Description

### DIGITAL GLASS SCALE - grey

Electronic bathroom scale with 6 mm thick tempered glass platform.  
Ultra slim design: only 19 mm height.  
4 high precision G sensors.  
Automatic on and off function.  
Lb/kg/st switchable by user.

LCD display 72x35 mm

Multilingual colour box and manual: GB, FR, IT, ES, DE, PL, PT, HU, RO, FI, SE, HR, GR, Arabic

## Technical Specifications

Capacity: max 150 kg (330 lb), min. 2 kg (4.4 lb)  
Sensitivity: 100 g (0.2 lb)  
Dimension: 300x300x19 mm  
Power: 1xCR2032 lithium battery - not included  
Warranty: 2 years