



Product Code	25725
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2449889
CND	V9003
EAN/UPC	8023279257250
UMDNS	14072
PARAF	903153031
GTIN	08023279257250
GMDN	67045

## Description

### GIMA TOURNIQUETS

- Model: STRIP-VELCRO Tourniquet
- Material: velcro
- Colour: green
- Length: 45 cm
- Width: 2.5 cm
- Fastening systems: Velcro+ring: use with one hand

Efficient and simple solutions for professional medical use.

Velcro "D" buckle: self applicable made of a nylon band threads through plastic "D" buckle

Multilingual box: GB, FR , IT, DE, ES , PT, GR, Arabic.