

DIABETIC BAG empty - nylon blue



Product Code	27255
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279272550
PARAF	901713976

Description

DIABETIC BAG - empty - blue

Thermal pharmaceutical bags designed specifically for carrying medicines that must be kept at cool temperatures, ideal for diabetics. Temperature pack included to provide a temperature of about 10°C for up to 10 hours. Exterior material is a water repellent vinyl nylon with a foam backing.

Size: 24x18xh 6 cm

Made in Italy.