



Product Code	27257
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	4012030006750

## Description

### SECA 807 DIGITAL SCALE

Digital scale with advanced technology for enhanced comfort.  
 Extra-flat design to easily step onto the platform.  
 Four specially constructed load cells made of stainless steel, guarantee high precision.  
 Supplied with power sensors, extralarge LCD numbers and high-quality glass surface.  
 Safety shock resistant glass construction.  
 Step-off, autohold and automatic switch ON/OFF functions.

Multilingual manual: GB, FR, IT, ES, PT, DE, PL, NL, GR, DK, SE, FI, NO, TR, JP, RU.

Digit height: 29 mm

## Technical Specifications

Capacity: 150 kg  
 Sensitivity: 100 g  
 Size: 310x310xh 37 mm  
 Weight: 2 kg  
 Batteries: 4xAAA 1,5V (not included)