

SPECULUM diam.3mm for 72560, W.A. series 20100, 20200, 20202, 21600, 21700,



Product Code	72570
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

EAN/UPC	0732094107975
UMDNS	13662
GMDN	34897

Description

WELCH-ALLYN KLEENSPEC DISPOSABLE EAR SPECULUM
Colour: black
Ø: 3 mm
To be used with: Welch-Allyn 20100, 20200, 20202, 21600, 21700
Minimum order: box of 100 (4 packs of 25 pcs.)