



Product Code	27258
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	5997307566801

Description

DIGITAL PET SCALE - medium

Precise full electronic multifunctional scale, designed for measuring your pet's weight available in two sizes.

Quality product with very competitive price.

Switchable kg, lb, st.

Minimal measuring weight: 3 kg.

Automatic switch-off.

Multilingual manual: GB, FR, IT, ES, PT, DE, PL, CZ, SE, FI, RO, GR, HR, BG.

Made in Europe.

Technical Specifications

Capacity/Sensitivity: 100 kg / 100 g

Display: LCD, 3"

Dimensions: 355x635 mm

Lithium battery: included, CR2032, 3 V