

# SECA 813 DIGITAL FLAT SCALE



|               |                   |
|---------------|-------------------|
| Product Code  | 27265             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 4012030020770 |

## Description

### SECA 813 DIGITAL FLAT SCALE

Electronic scale with high capacity up to 200 kg.  
Extra-wide platform and flat shape facilitate weighing of heavyweight patients and make stepping on very easy.  
Made of solid steel frame with tough rubber coating and non-slip rubber base, ensuring safety and stability.  
Kg/lbs/sts switch.  
Auto-hold, step-off and automatic switch ON/OFF functions.  
Multilingual manual: GB, FR, IT, ES, PT, DE, PL, NL, GR, DK, SE, FI, NO, TR, JP, RU.

## Technical Specifications

Capacity: 200 kg - 440 lbs  
Sensitivity: 100 g - 0,2 lbs  
Size: 433x373xh 47 mm  
Platform size: 400x300xh 45 mm  
Weight: 2,9 kg  
Batteries: 4xAAA 1,5V (not included)