

SOEHNLE 8320 FOLDABLE BABY SCALE



Product Code	27266
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	4006484830009
PARAF	987476557

Description

SOEHNLE 8320 FOLDABLE BABY SCALE

Foldable baby scale for easy transport.
Folding mechanism with magnetic lock.
Easy to open: remove the base and pull the handles apart.
Taring function: pads, blankets placed under the baby are tarable.
Steady weight display through filtering and hold function.
Large LCD display easy to read.
Switchable to kg/lb.

Manual in: DE, GB, FR, ES, IT, GR, CZ, Arabic.

Technical Specifications

Capacity: 15 kg
Sensitivity: 10 g
Dimensions: 395x315x105 mm (closed) - 740x315x53 mm (open)
Weight: 1.9 kg
Alkaline battery 2x1.5 V AAA - not included
Warranty: 2 years