

SECA 376 HOSPITAL DIGITAL BABY SCALE 20 kg Class III



Product Code	27267
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	76914
CND	Z12080402
EAN/UPC	4012030008181

Description

SECA 376 HOSPITAL SCALE - digital class III

- high-quality Seca scale technology with fine 5-gram graduation
- TARE, HOLD, auto-HOLD and BMIF for extra convenience
- extra large and ergonomic tray for safe weighing
- adjustable damping
- large display for easy reading

Technical Specifications

- Capacity: 20 kg
- Sensitivity: 5 g up to 7.5 kg, 10 g from 7.5 kg up
- Dimensions: 620 x 190 x 358 mm
- Weight: 3.7 kg
- Platform: 615 x 130 x 275 mm
- Battery operated