

SECA 354 DIGITAL BABY SCALE - 20 kg



Product Code	27282
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	4012030018586
PARAF	990604872

Description

SECA 354 - electronic

Two scales in one: with a comfortable tray to weigh in a reclining position, easily detachable to be converted into a floor scale. Breast-milk-intake function determines how much milk the baby has suckled.

Other functions: tare, auto-hold, kg/lbs switch over automatic switch off.

User manual: GB, FR, IT, ES, DE.

Technical Specifications

- Capacity: 20 kg
- Size: 552 x 156 x h 332 mm
- Battery operated
- Sensitivity: 10 g < 10 kg > 20 g
- Weight: 2.3 kg
- Made in Asia