



Product Code	27285
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No
RDM (NSIS)	1404831
CND	V0399
EAN/UPC	4012030012843

## Description

Traditional model with large easy-to-read dial. Pointer professional scale in strong metal structure white painted.

Consistent accuracy on any surface.  
Spacious platform design with no slip rubber mat, stable for safety.

Tare function.  
User manual in: GB, FR, IT, ES, DE, PT, PL.

## Technical Specifications

Capacity: 150 kg - 330 lbs  
Sensitivity: 500 g  $\pm$  2 g  
Size: 305x120x470 mm  
Weight: 3.5 kg