

SECA 799 DIGITAL SCALE with BMI Class III - 200 kg



Product Code	27292
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No
RDM (NSIS)	93791
CND	V0399
EAN/UPC	4012030007177
GMDN	109277

Description

SECA 799 - electronic with BMI - class III

Versatile, low maintenance and site independent with economic battery operation, ideal for daily use in hospital and doctors practices.

Column scale is easy to move around with its transport castors.

BMI function permits a reliable evaluation of the nutritional conditions of the patient.

Equipped with several fuctions: TARE, HOLD, weighing range switch-over, automatic switch off.

User manual in GB, FR, IT, ES, DE, PT, PL, NL, GR, DK, SE, FI, NO.

Technical Specifications

Capacity: 200 kg
Sensitivity: 100 g < 150 kg > 200 g
Weight: 6.3 kg
Dimensions: 294x831x417 mm
Platform: 272x280 mm